

Respireo **SOFT** Baby

Technical review


General presentation

Infant mask with vent system



Technical data

MINIMUM CONTINUOUS FLOW



Normative references	CE marking in accordance with European Directive 93/42/EEC, class IIa EN ISO 14971, ISO 5356-1, ISO 10993-1, EN ISO 17510-2
Resistance	at 50 lpm: 0,9 cmH2O - at 100 lpm: 3,7 cmH2O
Operating pressure	4-24 cm H2O
Temperature range	5°-40°C
Weight	54,5g (XS size with tube) - 31,7g (without tube)
Dimensions	86 mm x 70 mm x 298 mm
Weighted sound pressure A according to ISO3744	29,5 dB(A)
Dead space volume inside the mask	50 ml

Technical data

Codes used by supplier	KM225300 RESPIREO SOFT Baby - XS KM225400 RESPIREO SOFT Baby - S
Manufacturer	Air Liquide Medical Systems Spa
Mask type	Nasal mask
Patients indication	3 to 12 kgs (6-26 pounds)
Sizes of the mask	
Mask body (cushion)	1
Headgear	2 (XS, Small)
Forehead support	Yes (possibility to cut it, please refer to the process described in the UM)
Chin Support	No
Oxygen pressure port	Yes
Flexible mask body	Yes (entirely in silicone)
360 degree swivels	Yes
Vented version	Yes
Headgear adjustment points	5
Headgear attachment method	Hooks (slit)
Headgear included in the pack	Yes
Storage bag	Yes
User manual languages	Deutsch, English, Español, Français, Italiano, Greek, Netherlands, Swedish, Polski, Português, Chinese, Japanese, Russian, Czech
Requirements for humidifier	Optional
Made in	Italy
CE mark	Yes
Manufacturer's warranty	120 Days

Materials

COMPONENT	MATERIAL
Main body	Polycarbonate + Silicone
Tube assembly	Polycarbonate + Metalocen + Silicone +Acetalic
Headgear	Nylon + polyurethane